

BEAM

ToxiScan
Precision behind Purity



ToxiScan

Technical Specifications

Parameter	Value
Dimensions(L x B x H)	191 x 251 x 278 mm
Weight	7 Kg
Technology	FTIR
Sample Type	Liquid
Sample	Propylene Glycol and Glycerin
Adulterants	Diethylene glycol more than 1%
Sample Test Total Time	<5 mins
Sample Volume	14 μ L
Data Interface	USB
Display	4.3" Resistive Touchscreen
Warm up time	30 minutes
Power Supply	12V / 10A

☎ (+91) 9975568458

✉ info@beamoptics.in

➤ www.beamoptics.in

† PAP J188, Near Quality Circle Forum,
Telco Road, Bhosari MIDC, PCMC, Pune 411026.

BEAM
OPTICS